



# Anti-Methotrexate monoclonal antibody (DMAB-WZ0004)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	N/A
<b>Immunogen</b>	MTX-KLH
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	MTX04
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	N/A
<b>Size</b>	500 µg
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** METHOTREXATE; METHYLAMINOPTERIN; METHYLAMINOPTERINE; ANTIFOLAN; L-(+)-AMETHOPTERIN; 4-AMINO-N 10-METHYLPTEROYL-L-GLUTAMIC ACID