



Anti-Hep Par-1 monoclonal antibody, clone PDI2F6 (DMAB13904)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Hep Par-1 Monoclonal Antibody
Target	Hep Par-1
Immunogen	Formalin-fixed, rejected liver transplant.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	PDI2F6
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In cell culture supernant (BSA, sodium azide)
Preservative	See individual product datasheet
Storage	Store at 4°C.

BACKGROUND

Introduction The antigen recognized by the antibody has not been yet characterized. By native slot blot technique, this antibody localizes to the mitochondrial fraction of rat liver following sucrose

gradient centrifugation

Keywords

Hepatocyte Specific Antigen; Mouse Anti-Hep Par-1 Monoclonal Antibody; Anti-Hep Par-1 Monoclonal Antibody; Hep Par-1 Monoclonal Antibody Mouse Anti-Hep Par-1 MAb; Anti-Hep Par-1 MAb; Hep Par-1 MAb; Mouse Anti-Hep Par-1 Antibody; Anti-Hep Par-1 Antibody; H
