



Anti-Thyroglobulin Monoclonal antibody, clone N456 (DMAB31924)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-Thyroglobulin Monoclonal Antibody
Target	Thyroglobulin
Immunogen	Human Thyroglobulin
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	N456
Purification	Purified by ion exchange chromatography.
Conjugate	Unconjugated
Size	1 mg
Buffer	Supplied in PBS, pH 7.4 with 0.05% sodium azide (NaN ₃) added as preservative.
Preservative	0.05% Sodium Azide
Storage	Store antibody at 2-8 °C. Aliquot and freeze at -20 °C for long term storage.

BACKGROUND

Introduction	Precursor of the iodinated thyroid hormones thyroxine (T ₄) and triiodothyronine (T ₃).
---------------------	---

Keywords

AITD4antibody; AITD3 antibody; hTG antibody; TDH3 antibody; Tg antibody; Tgn antibody;
THYG_HUMAN antibody; Thyroglobulin antibody; Mouse Anti-Thyroglobulin Monoclonal
Antibody; Anti-Thyroglobulin Monoclonal Antibody; Thyroglobulin Monoclonal Antibody; Mo

GENE INFORMATION

Entrez Gene ID[7038](#)**UniProt ID**[P01266](#)