



Anti-Norovirus Capsid protein Monoclonal antibody, Clone C115M (DMAB4017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes the capsid protein of Norovirus (native and recombinant). Does not cross-react with Rotavirus, Adenovirus or Astrovirus.
Target	Norovirus Capsid protein
Immunogen	Purified native Norwalk virus, strain 8FIIa
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Norovirus
Clone	C115M
Affinity Constant	Not determined
Purification	90% pure. Protein A chromatography
Conjugate	Unconjugated
Applications	<p>Suitable for use in ELISA for detection in stool. Pairs with Catalog #C01257M in sandwich immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p> <p>Recommended pairs for sandwich immunoassay:</p> <ul style="list-style-type: none"> • Capture DMAB4017 DMAB4016

- **Detection**
[DMAB4016](#)
[DMAB4017](#)

Suggested pair for testing (Capture - Detection): DMAB4017 - [DMAB4016](#)

Format	Purified, Liquid
Concentration	100ug/ml (OD280nm, E0.1% = 1.3)
Size	1 mg
Buffer	0.01M PBS, pH 7.2 This product contains no stabilizing proteins
Preservative	0.1% Sodium Azide
Storage	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Norovirus is an RNA virus of the Caliciviridae taxonomic family, causes approximately 90% of epidemic nonbacterial outbreaks of gastroenteritis around the world and is responsible for 50% of all foodborne outbreaks of gastroenteritis in the US. Norovirus affects people of all ages. The viruses are transmitted by faecally contaminated food or water and by person to person contact.
Keywords	Caliciviridae; Norovirus; VP1; Group IV