



Anti-PIV type 3 Hemagglutinin Monoclonal antibody, Clone C3292M (DMAB4035)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific for hemagglutinin of Parainfluenza virus, type 3. Also recognizes Bovine Parainfluenza virus, type 3.
Target	PIV type 3 Hemagglutinin
Immunogen	Parainfluenza Type 3 viral lysate
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	PIV
Clone	C3292M
Affinity Constant	Not determined
Purification	90% pure. Protein A chromatography
Conjugate	Unconjugated
Applications	<p>Suitable for use in ELISA and IFA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p> <p>Recommended pairs for sandwich immunoassay:</p> <ul style="list-style-type: none"> • Capture DMAB4035 DMAB4036 • Detection

[DMAB4036](#)

[DMAB4035](#)

Suggested pair for testing (Capture - Detection): DMAB4035 - [DMAB4036](#)

Format	Purified, Liquid
Concentration	100ug/ml (OD280nm, E0.1% = 1.3)
Size	1 mg
Buffer	PBS, pH 7.2. This product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Storage	Upon receipt, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Parainfluenzahemagglutinin-neuraminidase is a type of Hemagglutinin-neuraminidase produced by Parainfluenza
Keywords	HPIV3; Parainfluenza 3 virus; PIV3; Parainfluenza 3, Hemagglutinin