



Mouse Anti-Human IL6 Monoclonal Antibody, clone 20 (DMABB-JX195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant human IL-6 protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	20
Purification	Protein A/G Purified
Conjugate	Unconjugated
Applications	IHC, WB Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	50 µl, 100 µl, 200 µl
Buffer	PBS
Preservative	None
Storage	Store at -80°C. Avoid freeze/thaw cycles.

Ship

Dry ice

BACKGROUND

Introduction

Interleukin 6 (IL-6) is an interleukin that acts as both a pro-inflammatory cytokine and an antiinflammatory myokine. In humans, it is encoded by the IL6 gene.

Keywords

IL6; Interleukin-6; BSF2; HSF; IFNB2

GENE INFORMATION

Official Symbol

IL-6

UniProt ID

[P05231](#)
