



Rabbit Anti-Human IL-6 monoclonal antibody, clone S370U0 (DMABB-JX365)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant IL-6 expressed in mammalian cells
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	S370U0
Purification	Protein A/G purified
Conjugate	Unconjugated
Applications	ELISA, LFIA, CLIA We recommend the following as antibody pair for CLIA (Capture - Detection): DMABB-JX362 - DMABB-JX365 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None

Storage

Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

IL-6 is a potent lymphoid cell growth factor that stimulates the growth and survival of certain B cells and T cells. IL-6 plays a role in host defense, acute phase reactions, immune response, and hematopoiesis. IL-6 is expressed by T cells, B cells, monocytes, fibroblasts, hepatocytes, endothelial cells, and keratinocytes.

Keywords

IL6;interleukin 6;HGF;HSF;BSF2;IL-6;IFNB2;interleukin-6;CDF;BSF-2;IFN-beta-2;interferon beta-2;interleukin BSF-2;interferon;beta 2;hybridoma growth factor;CTL differentiation factor;Bcell stimulatory factor 2;B-cell differentiation factor;PF-04236921

GENE INFORMATION

Official Symbol

IL-6

Synonyms

IL6; interleukin 6; HGF; HSF; BSF2; IL-6; IFNB2; interleukin-6; CDF; BSF-2; IFN-beta-2; interferon beta-2; interleukin BSF-2; interferon; beta 2; hybridoma growth factor; CTL differentiation factor; B-cell stimulatory factor 2; B-cell differentiation factor; PF-04236921

Entrez Gene ID

[3569](#)

UniProt ID

[P05231](#)