



Mouse Anti-Human h-FABP monoclonal antibody, clone 38D2D (DMABB-JX580)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	h-FABP
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	38D2D
Purification	Protein A/G purified
Conjugate	Unconjugated
Applications	ELISA, CLIA We recommend the following as antibody pair (Capture - Detection): DMABB-JX580 - DMABB-JX581 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4, containing 0.03% Proclin 300
Preservative	0.03% Proclin 300
Storage	Store at -20°C or below. Avoid freeze / thaw cycles.

BACKGROUND

Keywords

FABP3; fatty acid binding protein 3, muscle and heart; MDGI; FABP11; H-FABP; M-FABP; O-FABP; fatty acid-binding protein, heart; fatty acid binding protein 11; mammary-derived growth inhibitor; muscle fatty acid-binding protein; Fatty acid-binding protein
