



## Human Anti-Human h-FABP monoclonal antibody, clone ID30 (DMABB-JX581)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	h-FABP
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	ID30
<b>Purification</b>	Protein A/G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, CLIA We recommend the following as antibody pair (Capture - Detection): DMABB-JX580 - DMABB-JX581 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS, pH 7.4, containing 0.03% Proclin 300
<b>Preservative</b>	0.03% Proclin 300
<b>Storage</b>	Store at -20°C or below. Avoid freeze / thaw cycles.

# BACKGROUND

**Keywords**

FABP3; fatty acid binding protein 3, muscle and heart; MDGI; FABP11; H-FABP; M-FABP; O-FABP; fatty acid-binding protein, heart; fatty acid binding protein 11; mammary-derived growth inhibitor; muscle fatty acid-binding protein; Fatty acid-binding protein

---