



# Human Anti-Sebastes marinus Seb m 1 monoclonal antibody, clone B83G (DMABC-JX164)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	Sebastes marinus Seb m 1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Sebastes marinus
<b>Clone</b>	B83G
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IA, Control. This recombinant human antibody can be used as an alternative to serum based calibrators and controls. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid, Purified
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µg
<b>Buffer</b>	PBS or Tris-Gly
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Ship

Wet ice

---

## BACKGROUND

Keywords

Sebastes marinus; Seb m 1

---