



Human Anti-Rana esculenta Ran e 2 monoclonal antibody, clone B64G (DMABC-JX187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Rana esculenta Ran e 2
Isotype	IgG
Source/Host	Human
Species Reactivity	Rana esculenta
Clone	B64G
Conjugate	Unconjugated
Applications	IA, Control. This recombinant human antibody can be used as an alternative to serum based calibrators and controls. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid, Purified
Concentration	Lot specific
Size	100 µg
Buffer	PBS or Tris-Gly
Preservative	None
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Ship

Wet ice

BACKGROUND

Keywords

Rana esculenta; Ran e 2
