



## Rabbit Anti-Human SLC38A2 monoclonal antibody, clone D7 (DMABC-JX273)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	SNAT2
<b>Immunogen</b>	SNAT2 fusion protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Clone</b>	D7
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB 1:5000-1:50000; IHC 1:4000-1:16000; ELISA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	50 µl, 100 µl
<b>Buffer</b>	PBS with 0.02% sodium azide and 50% glycerol
<b>Preservative</b>	0.02% sodium azide

**Storage**

Store at -20°C.

---

## **GENE INFORMATION**

**Entrez Gene ID**

[54407](#)

---

**UniProt ID**

[Q96QD8](#)

---