



## Mouse Anti-S. abortusequi ompA monoclonal antibody, clone 3A6A2 (DMABC-JX293)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Target</b>	S. abortusequi ompA
<b>Immunogen</b>	Recombinant protein ompA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Salmonella Abortusequi
<b>Clone</b>	3A6A2
<b>Purification</b>	Affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 ea
<b>Buffer</b>	PBS
<b>Preservative</b>	< 0.1% sodium azide

**Storage**

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

---

## BACKGROUND

**Keywords**

Salmonella abortusequi; S. abortusequi; Salmonella abortus equi; S. abortus equi; S. abortusequi ompA; S. abortus equi ompA; ompA

---