



Anti-HBV Envelope Antigen Monoclonal Antibody [HRP], clone N21182033 (DMABH008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Matched pair antibodies available for HBeAg antibody: DMABH009
Specificity	HBV
Immunogen	HBeAg antibody (HRP) was raised in mouse using hepatitis B virus as the immunogen.
Source/Host	Mouse
Species Reactivity	To be determined by end-user
Clone	N21182033
Conjugate	HRP
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS and 50% glycerol, pH 7.4.
Preservative	None
Storage	Store at 4°C for short term storage. Aliquot and store at -20°C for long term storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Hepatitis B Virus Envelope Antigen; hepatitis B e antigen; HBeAg; HBV; Hepatitis B Virus; Hepatitis B 'E' Antigen; HBV Envelope Antigen