



Anti-HBV Envelope Antigen Monoclonal Antibody, clone N1182032 (DMABH009)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DMABH009 can be used as a detection antibody in a matched pair immunoassay in conjunction with DMABH008 as the capture antibody
Specificity	HBV
Immunogen	HBeAg antibody was raised in mouse using hepatitis B Virus as the immunogen.
Source/Host	Mouse
Species Reactivity	To be determined by end-user
Clone	N1182032
Conjugate	Unconjugated
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH 7.4.
Preservative	None
Storage	Store at 4°C for short term storage. Aliquot and store at -20°C for long term storage. Avoid repeated freeze/thaw cycles.

BACKGROUND

Keywords Hepatitis B Virus Envelope Antigen; hepatitis B e antigen; HBeAg; HBV; Hepatitis B Virus; Hepatitis B 'E' Antigen; HBV Envelope Antigen