



Magic™ Anti-Cancer Antigen 15-3 (CA 15-3) MUC 1 Monoclonal Antibody, clone IN265 (DMABP-L02)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	MUC1
Immunogen	Purified native MUC1
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	IN265
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Format	Ion exchange purified monoclonal antibody
Size	1 mg
Buffer	10 mM Tris, pH 7.5 with 0.15M NaCl
Preservative	0.05% Sodium Azide
Storage	2-8°C

GENE INFORMATION

Entrez Gene ID [4582](#)

UniProt ID [P15941](#)
