



Magic™ Anti-Canine Heartworm Monoclonal Antibody, clone IN274 (DMABP-L04)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Dirofilaria immitis
Target	Canine Heartworm
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Dog
Clone	IN274
Conjugate	Unconjugated
Applications	ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	2-8°C