



# Magic™ Anti-Dengue Virus (NS1) Monoclonal Antibody, clone IN275 (DMABP-L07)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	NS1 protein of all 4 serotypes
<b>Target</b>	DENV
<b>Immunogen</b>	Recombinant Dengue NS1
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	N/A
<b>Clone</b>	IN275
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, IFA, LFIA
<b>Format</b>	Purified Monoclonal Antibody
<b>Size</b>	1 mg
<b>Buffer</b>	0.01 M Phosphate Buffered Saline, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	2-8°C