



Magic™ Anti-Dengue Virus Monoclonal Antibody, clone IN037 (DMABP-L08)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	NS1 protein of all 4 serotypes
Target	DENV
Immunogen	Recombinant Dengue NS1
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N/A
Clone	IN037
Conjugate	Unconjugated
Applications	LFIA, ELISA, IFA
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Avoid repeated freezing and thawing.