



# Magic™ Anti-Distemper Virus, Canine Monoclonal Antibody, clone IN207 (DMABP-L10)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	CDV surface envelope antigen
Target	CDV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Dog
Clone	IN207
Conjugate	Unconjugated
Applications	LFIA, ELISA
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	2-8°C