



Magic™ Anti-Feline Leukemia Virus (p27) Monoclonal Antibody, clone IN200 (DMABP-L12)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	FeLV p27
Target	FeLV P27
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Cat
Clone	IN200
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Protein A Purified Monoclonal Antibody
Size	500 µg
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid multiple freeze-thaw cycles