



Magic™ Anti-Human Immunoglobulin E (IgE) Monoclonal Antibody, clone GB429 (DMABP- L14)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	IgE
Immunogen	Human IgE
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	GB429
Conjugate	Unconjugated
Applications	ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.05 M Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	2-8°C