



# Magic™ Anti-Human Immunoglobulin E Monoclonal Antibody, clone IN228 [HRP] (DMABP-L15)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	IgE
Immunogen	Human IgE
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	IN228
Conjugate	HRP
Applications	ELISA
Format	Purified Monoclonal Antibody-HRP Conjugate
Size	1 mg
Buffer	Protein stabilization buffer, pH 7.6
Preservative	None
Storage	4-8°C