



# Magic™ Anti-Listeria monocytogenes Monoclonal Antibody, clone IN328 (DMABP- L17)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                          |
|--------------------|------------------------------------------|
| Specificity        | L. monocytogenes p60                     |
| Target             | L. monocytogenes                         |
| Isotype            | IgG2b                                    |
| Source/Host        | Mouse                                    |
| Species Reactivity | L. monocytogenes                         |
| Clone              | IN328                                    |
| Conjugate          | Unconjugated                             |
| Applications       | ELISA, IFA                               |
| Format             | Purified Monoclonal Antibody             |
| Size               | 1 mg                                     |
| Buffer             | 0.01 M Phosphate Buffered Saline, pH 7.2 |
| Preservative       | 0.1% Sodium Azide                        |
| Storage            | Short Term: 2-8°C Long Term: -20°C       |