



# Magic™ Anti-Malaria (HRP-2) Monoclonal Antibody, clone IN098 (DMABP-L20)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	P. falciparum LDH
<b>Immunogen</b>	Recombinant Malaria HRP-2 antigen
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	P. falciparum
<b>Clone</b>	IN098
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	LFIA, ELISA
<b>Format</b>	Purified Monoclonal Antibody
<b>Size</b>	1 mg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Short Term ( $\leq 2$ weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles.