



Magic™ Anti-Malaria (pLDH) Monoclonal Antibody, clone IN090 (DMABP-L23)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	P. falciparum LDH
Immunogen	Recombinant Malaria pLDH antigen
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	P. falciparum
Clone	IN090
Conjugate	Unconjugated
Applications	LFIA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	None
Storage	Short Term (≤ 2 weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles.