



# Magic™ Anti-Myoglobin (Mb) Monoclonal Antibody, clone IN323 (DMABP-L25)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	MB
Immunogen	Native Human Myoglobin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN323
Conjugate	Unconjugated
Applications	LFIA
Format	Liquid; Purified Monoclonal Antibody
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	None
Storage	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing

## GENE INFORMATION

Entrez Gene ID [4151](#)

UniProt ID

[P02144](#)

---