



Magic™ Anti-Prolactin Monoclonal Antibody, clone GB420 (DMABP-L40)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Prolactin
Immunogen	Human Prolactin
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	GB420
Conjugate	Unconjugated
Applications	ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.05 M Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	Short Term (≤ 2 weeks): Refrigerated at 2-8°C. Long Term: Frozen at -20°C. Avoid repeated freeze/thaw cycles.

GENE INFORMATION

Entrez Gene ID [5617](#)

UniProt ID

[P01236](#)
