



# Magic™ Anti-Prolactin Monoclonal Antibody, clone IN225 [HRP] (DMABP-L42)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Prolactin
Immunogen	Human Prolactin
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	IN225
Conjugate	HRP
Applications	ELISA
Format	Purified Monoclonal Antibody-HRP Conjugate
Size	1 mg
Buffer	protein stabilization buffer, pH 7.6
Preservative	None
Storage	4-8°C

## GENE INFORMATION

Entrez Gene ID [5617](#)

UniProt ID

[P01236](#)

---