



Magic™ Anti-Respiratory Syncytial Virus (RSV) Monoclonal Antibody, clone K469 (DMABP-L43)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Fusion protein of RSV types A & B
Target	RSV
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	N/A
Clone	K469
Conjugate	Unconjugated
Applications	LFIA, ELISA, IFA
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Avoid repeated freezing and thawing.