



# Magic™ Anti-Salmonella LPS Core Monoclonal Antibody, clone IN050 (DMABP-L45)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	LPS core
Target	Salmonella
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Salmonella
Clone	IN050
Conjugate	Unconjugated
Applications	LFIA, ELISA, IFA
Format	Protein A Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Avoid repeated freezing and thawing.