



Magic™ Anti-Streptococcus pneumoniae Monoclonal Antibody, clone IN299 (DMABP- L47)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	S. pneumoniae common cell wall polysaccharide (urinary antigen)
Target	S. pneumoniae
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Clone	IN299
Conjugate	Unconjugated
Applications	LFIA, ELISA, IFA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	0.01 M Phosphate Buffered Saline, pH 7.2
Preservative	0.1% Sodium Azide
Storage	2-8°C