



Magic™ Anti-Thyroid Stimulating Hormone (TSH) Monoclonal Antibody, clone GB406 (DMABP-L48)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	TSH
Immunogen	Human TSH β -subunit
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	GB406
Conjugate	Unconjugated
Applications	ELISA
Format	Purified Monoclonal Antibody
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.1% Sodium Azide
Storage	2-8°C

GENE INFORMATION

Entrez Gene ID [1081](#)

UniProt ID [P01215](#)
