



Anti-GMEB1 monoclonal antibody, clone 3B9 (DMABT-H13834)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-GMEB1 Monoclonal Antibody
Target	GMEB1
Immunogen	GMEB1 (NP_077808, 467 a.a. ~ 563 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3B9
Conjugate	Unconjugated
Applications	WB, IHC, IF, sELISA, ELISA
Sequence Similarities	TAMQDGSTLGNMTTMVSPVELVAMESGLTSAIQAVESTSEDGQTIIIDPAPDPEAEDTE GKAVILETELRTTEKVVAAEMEEHQHQVHNVEIVVLED
Buffer	In 1x PBS, pH 7.2
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GMEB1 glucocorticoid modulatory element binding protein2[Homo sapiens]
Official Symbol	GMEB1
Synonyms	GMEB1; glucocorticoid modulatory element binding protein 1; glucocorticoid modulatory element-binding protein 1; P96PIF; PIF96; GMEB-1; PIF p96; DNA-binding protein p96PIF; parvovirus initiation factor p96;
Entrez Gene ID	10691
Protein Refseq	NP_006573
UniProt ID	B2RCN8
Chromosome Location	1p35
Function	DNA binding; binding; metal ion binding; sequence-specific DNA binding transcription factor activity; transcription coactivator activity;