



Anti-EOMES monoclonal antibody, clone 32Nbht9 (DMABT-H22525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-EOMES Monoclonal Antibody
Target	EOMES
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Human
Clone	32Nbht9
Conjugate	Unconjugated
Applications	FC, WB
Size	25 µg
Preservative	None
Storage	-20°C, avoid repeated freeze/thaw cycles

GENE INFORMATION

Gene Name	EOMES eomesodermin [Homo sapiens]
Official Symbol	EOMES
Synonyms	EOMES; eomesodermin; eomesodermin (Xenopus laevis) homolog; eomesodermin homolog; T box brain2; TBR2; TBR-2; T-brain-2; T-box brain2; t box, brain, 2; T-box brain protein 2;

Entrez Gene ID	8320
Protein Refseq	NP_005433
UniProt ID	B7Z4K0
Chromosome Location	3p24.1
Pathway	Downstream signaling in naive CD8+ T cells, organism-specific biosystem; IL12-mediated signaling events, organism-specific biosystem.
Function	DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;