



## Anti-*S. marcescens* Monoclonal antibody, Clone D/P6P (DMABT-Z59999)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	Mouse Anti-Serratia marcescens Monoclonal Antibody
<b>Target</b>	<i>S. marcescens</i>
<b>Immunogen</b>	Native Serratia marcescens (ATCC 13880)
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	<i>S. marcescens</i>
<b>Clone</b>	D/P6P
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, ICC/IF
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	PBS
<b>Preservative</b>	10mM Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

# BACKGROUND

## Introduction

*Serratia marcescens* is a species of Gram negative bacteria in the family Enterobacteriaceae. It is a broad host range pathogen, and is also capable of opportunistic infections of humans. Some strains of *S. marcescens* are capable of producing a pigment cal

---

## Keywords

*Serratia marcescens*

---