



Mouse Anti-EBV EBNA2 Monoclonal antibody, Clone RG4 (DMABT-Z60760)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-EBNA2 Monoclonal Antibody
Target	EBV
Immunogen	Recombinant fusion protein against EBNA2.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	EBV
Clone	RG4
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, IP, IHC-Fr, ICC, IF
Cellular Localization	Nuclear matrix of host cells
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The Epstein-Barr Virus (EBV), also called Human herpesvirus5(HHV-4), is a virus of the herpes family and is one of the most common viruses in humans. EBNA2 is one of the few genes of Epstein-Barr virus which are necessary for immortalization of human prim
Keywords	EBNA2; EBV; EBV nuclear antigen 2; Epstein Barr nuclear antigen 2; HHV4; Human Herpesvirus 4;