



Anti-Tetanus Toxoid polyclonal antibody (DMABT007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti-Tetanus Toxoid polyclonal antibody.
Specificity	Reacts with Tetanus Toxoid.
Immunogen	Immunogen was purified Tetanus Toxoid.
Source/Host	Goat
Species Reactivity	Bacteria
Clone	D46H48
Purification	The IgG fraction is purified from goat anti-serum to >95% purity.
Conjugate	Unconjugated
Applications	ELISA, WB, IF, IP
Size	1 ml
Buffer	Supplied in PBS, pH 7.2, with 0.1% sodium azide added as preservative.
Preservative	0.1% Sodium Azide
Storage	Store antibody at -20°C. Avoid freeze-thaw cycles.

BACKGROUND

Introduction	Tetanus toxin is an extremely potent neurotoxin produced by the vegetative cell of Clostridium tetani in anaerobic conditions, causing tetanus. It has no known function for clostridia in the soil
--------------	---

environment where they are normally encountered. It is also

Keywords

Clostridium tetani toxoid; Tentoxylisin; Tetanus toxin; tetX; TETANUS TOXOID,
CLOSTRIDIUM TETANI
