



# Anti-C. albicans polyclonal antibody (DMABT009)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Polyclonal to Candida Albicans antibody. The antibody (Cat# DMABT009) is suitable for ELISA, Western blot, Immunofluorescence, Immunohistochemistry.
<b>Specificity</b>	Reacts with multiple Candida albicans proteins. Some cross reactivity with other yeast.
<b>Immunogen</b>	Immunogen was prepared from a whole cell lysate of Candida albicans.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Purification</b>	The IgG fraction is purified from rabbit anti-serum to >95% purity.
<b>Conjugate</b>	Unconjugated
<b>Size</b>	1 ml
<b>Buffer</b>	Supplied in PBS, pH 7.2, with 0.1% sodium azide added as preservative.
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store antibody at -20°C. Avoid freeze-thaw cycles.