



Anti-VASH1 (aa 3-100) polyclonal antibody (DPAB-DC1005)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	VASH1 (vasohibin 1) is a protein-coding gene. Diseases associated with VASH1 include female breast carcinoma, and severe pre-eclampsia. An important paralog of this gene is VASH2.
Immunogen	KIAA1036 (NP_055724, 3 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is GGKKVAGGGSSGATPTSAATAPSGVRRLETSEG TSAQRDEEPEEEGEEDLRDGGVPFFV NRGGLPVDEATWERMWKHVAKIHPDGEKVAQRIRGATD
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	VASH1 vasohibin 1 [Homo sapiens (human)]
Official Symbol	VASH1

Synonyms	VASH1; vasohibin 1; KIAA1036; vasohibin-1;
Entrez Gene ID	22846
Protein Refseq	NP_055724
UniProt ID	Q7L8A9
Chromosome Location	14q24.3