



Anti-BAHD1 polyclonal antibody (DPAB-DC1007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BAHD1 (bromo adjacent homology domain containing 1) is a protein-coding gene. Diseases associated with BAHD1 include listeriosis. GO annotations related to this gene include chromatin binding and DNA binding.
Specificity	Specific to recombinant protein GX0582. This antibody detects mBAHD1 protein.
Immunogen	Recombinant GST fusion protein corresponding to 170 mouse Bahd1. The sequence is AIRKSYQAVERRHGETIRVRDVLKSGPRKTSTPYVAKISALWENPESGELMMSLLWYYR PEHLQGGRSPSMHEPLQNEVFASRHQDQNSVACIEEKCYVLTFAEYCRFCAMAKRRGEGL PSRKTALVPPSADYSTPPHRTVPEDTDPELVFLCRHVYDFRHGRILKNPQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Bahd1 bromo adjacent homology domain containing 1 [Mus musculus (house mouse)]
Official Symbol	BAHD1
Synonyms	BAHD1; bromo adjacent homology domain containing 1; Gm117; AL022997; AW541238; mKIAA0945; bromo adjacent homology domain-containing 1 protein; BAH domain-containing protein 1;
Entrez Gene ID	228536
Protein Refseq	NP_001038988
UniProt ID	Z4YJL0
Chromosome Location	2 E5; 2
Pathway	XPodNet - protein-protein interactions in the podocyte expanded by STRING;
Function	chromatin binding;