



Anti-ELL2 (aa 253-352) polyclonal antibody (DPAB-DC1017)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ELL2 (elongation factor, RNA polymerase II, 2) is a protein-coding gene. Diseases associated with ELL2 include eaf, and hepatocellular carcinoma. An important paralog of this gene is MARVELD2.
Immunogen	ELL2 (NP_036213, 253 a.a. ~ 352 a.a) partial recombinant protein with GST tag. The sequence is SKDLSYTLKDYVFKELQRDWPGYSEIDRRSLESVLSRKLNPSQATGTSRSESPVCSSRD AVSSPQKRLLDSEFIDPLMNKKARISHLTNRVPPTLNHGL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ELL2 elongation factor, RNA polymerase II, 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	ELL2
Synonyms	ELL2; elongation factor, RNA polymerase II, 2; RNA polymerase II elongation factor ELL2; ELL-related RNA polymerase II, elongation factor;
Entrez Gene ID	22936
Protein Refseq	NP_036213
UniProt ID	O00472
Chromosome Location	5q15