



Anti-ZCCHC11 polyclonal antibody (DPAB-DC1038)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZCCHC11 (zinc finger, CCHC domain containing 11) is a protein-coding gene. GO annotations related to this gene include RNA uridylyltransferase activity and nucleic acid binding. An important paralog of this gene is ENSG00000255508.
Specificity	Specific to recombinant protein GX2516. This antibody detects mZCCHC11 protein.
Immunogen	Recombinant GST fusion protein corresponding to 249 amino acids of mouse Zcchc11.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Zcchc11 zinc finger, CCHC domain containing 11 [Mus musculus (house mouse)]
Official Symbol	ZCCHC11
Synonyms	ZCCHC11; zinc finger, CCHC domain containing 11; PPAPD3; mKIAA0191; 6030404K05Rik; 9230115F04Rik; terminal uridylyltransferase 4; TUTase 4; zinc finger CCHC domain-containing protein 11;
Entrez Gene ID	230594
Protein Refseq	NP_780681
UniProt ID	B2RX14
Chromosome Location	4 C7; 4
Function	RNA uridylyltransferase activity; metal ion binding; nucleic acid binding; nucleotidyltransferase activity