



Anti-TRIM35 (aa 394-493) polyclonal antibody (DPAB-DC1043)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified.
Immunogen	TRIM35 (NP_055881, 394 a.a. ~ 493 a.a) partial recombinant protein with GST tag. The sequence is SGFWYVCRTQGVGDHCVTSDPATSPVLVAIPRRRLRVELECEEGELSFYDAERHCHLYTF HARFGEVRPYFYLGARGAGPPEPLRICPLHISVKEELD
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [TRIM35 tripartite motif containing 35 \[Homo sapiens \(human\) \]](#)

Official Symbol	TRIM35
Synonyms	TRIM35; tripartite motif containing 35; HLS5; MAIR; tripartite motif-containing protein 35; tripartite motif-containing 35; hemopoietic lineage switch protein 5;
Entrez Gene ID	23087
Protein Refseq	NP_741983
UniProt ID	Q9UPQ4
Chromosome Location	8p21.2
Function	molecular_function; zinc ion binding;