



Anti-FAM208A polyclonal antibody (DPAB-DC1071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM208A (family with sequence similarity 208, member A) is a protein-coding gene. Diseases associated with FAM208A include retinoblastoma, and cutaneous t cell lymphoma. An important paralog of this gene is FAM208B.
Immunogen	Recombinant protein corresponding to amino acids of recombinant C3orf63. The sequence is TIIKASAKGGNLPPVSPND SGAKIASNPLERHVIPVSSSDFNKHLLEPLCSDPLKDTN SDEQHSTSALTEVEMNQPQHATELMVTSDHIVPGDMAREPVEETTKSPSDVNISAQPALS N
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and

thawing.

GENE INFORMATION

Gene Name	FAM208A family with sequence similarity 208, member A [Homo sapiens (human)]
Official Symbol	FAM208A
Synonyms	FAM208A; family with sequence similarity 208, member A; RAP140; se89-1; C3orf63; protein FAM208A; CTCL tumor antigen se89-1; retinoblastoma-associated protein 140; retinoblastoma-associated protein RAP140;
Entrez Gene ID	23272
Protein Refseq	NP_001106207
UniProt ID	Q9UK61
Chromosome Location	3p14.3
Function	poly(A) RNA binding;