



Anti-SASH1 (aa 1066-1175) polyclonal antibody (DPAB-DC1078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SASH1 (SAM and SH3 domain containing 1) is a protein-coding gene. Diseases associated with SASH1 include cleft lip, and colon cancer. An important paralog of this gene is SASH3.
Immunogen	SASH1 (NP_056093, 1066 a.a. ~ 1175 a.a) partial recombinant protein with GST tag. The sequence is ENTSLQEHGVKLGPAITRKVSCARGVDLETLTENKLHAEGIDLTEEPYSDKHGRCGIPEA LVQRYAEDLDQPERDVAANMDQIRVKQLRKQHRMAIPSGGLTEICRKPVS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SASH1 SAM and SH3 domain containing 1 [Homo sapiens (human)]
Official Symbol	SASH1

Synonyms	SASH1; SAM and SH3 domain containing 1; SH3D6A; dJ323M4; dJ323M4.1; SAM and SH3 domain-containing protein 1; 2500002E12Rik; proline-glutamate repeat-containing protein;
Entrez Gene ID	23328
Protein Refseq	NP_056093
UniProt ID	O94885
Chromosome Location	6q24.3
Function	mitogen-activated protein kinase kinase kinase binding; protein C-terminus binding; protein complex scaffold; protein kinase binding