



Anti-VPS8 (aa 1329-1428) polyclonal antibody (DPAB-DC1084)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	VPS8 (vacuolar protein sorting 8 homolog (<i>S. cerevisiae</i>)) is a protein-coding gene. GO annotations related to this gene include zinc ion binding.
Immunogen	KIAA0804 (NP_001009921, 1329 a.a. ~ 1428 a.a) partial recombinant protein with GST tag. The sequence is RITPSQVKMSPSYHQSKGDPTAKKGTSEPVLDPPQIQAFDQLCRLYRGSSRLALLTELSQ NRSSSEYRPFSGSQSAPAFNSIFQNFQLQLIPPPVTED
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	VPS8 vacuolar protein sorting 8 homolog (<i>S. cerevisiae</i>) [<i>Homo sapiens</i> (human)]
Official Symbol	VPS8

Synonyms	VPS8; vacuolar protein sorting 8 homolog (S. cerevisiae); KIAA0804; vacuolar protein sorting-associated protein 8 homolog;
Entrez Gene ID	23355
Protein Refseq	NP_001009921
UniProt ID	B3KPR6
Chromosome Location	3q27.2
Function	zinc ion binding;