



Anti-ZNF281 (aa 787-893) polyclonal antibody (DPAB-DC1115)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNF281 (zinc finger protein 281) is a protein-coding gene. GO annotations related to this gene include transcription corepressor activity and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is ZNF148.
Immunogen	ZNF281 (NP_036614, 787 a.a. ~ 893 a.a) partial recombinant protein with GST tag. The sequence is KLTSQKEQKNLESSTGFQIPSQELASQIDPQKDIEPRTTYQIENFAQAFGSQFKSGSRVP MTFITNSNGEVDHRVRTSVSDFSGYTNMMSDVSEPCSTRVKTPTSQS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ZNF281 zinc finger protein 281 \[Homo sapiens \(human\) \]](#)

Official Symbol	ZNF281
Synonyms	ZNF281; zinc finger protein 281; ZBP-99; ZNP-99; ZNP-99 transcription factor; transcription factor ZBP-99; zinc finger DNA-binding protein 99; GC-box-binding zinc finger protein 1;
Entrez Gene ID	23528
Protein Refseq	NP_001268222
UniProt ID	Q9Y2X9
Chromosome Location	1q32.1
Function	core promoter binding; metal ion binding; protein binding; sequence-specific DNA binding