



Anti-FRRS1L (aa 158-267) polyclonal antibody (DPAB-DC1135)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FRRS1L (ferric-chelate reductase 1-like) is a protein-coding gene, and is affiliated with the lncRNA class. Diseases associated with FRRS1L include dysautonomia. GO annotations related to this gene include molecular_function. An important paralog of this gene is FRRS1.
Immunogen	C9orf4 (NP_055149, 158 a.a. ~ 267 a.a) partial recombinant protein with GST tag. The sequence is FRYGKPGCNAETCDYFLSYRMIGADVEFELSADTDGWVAVGFSSDKKMGDDVMACVHDD NGRVRIQHFYNVGQWAKEIQRNPARDEEGVFENNRVTCRFKRPVNVPRDE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [FRRS1L ferric-chelate reductase 1-like \[Homo sapiens \(human\) \]](#)

Official Symbol	FRRS1L
Synonyms	FRRS1L; ferric-chelate reductase 1-like; CG6; CG-6; C9orf4; DOMON domain-containing protein FRRS1L; brain protein CG-6; ferric-chelate reductase 1-like protein;
Entrez Gene ID	23732
Protein Refseq	NP_055149
UniProt ID	Q9P0K9
Chromosome Location	9q31
Function	molecular_function;