



Anti-TSPYL5 polyclonal antibody (DPAB-DC1142)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TSPYL5 (TSPY-like 5) is a protein-coding gene. Diseases associated with TSPYL5 include lung squamous cell carcinoma, and acute myeloid leukemia. GO annotations related to this gene include protein binding. An important paralog of this gene is SET.
Specificity	Specific to recombinant protein GX1663. This antibody detects mTSPYL5 protein. It also recognizes human TSPYL5 protein.
Immunogen	Recombinant GST fusion protein corresponding to 245 amino acids of mouse Tspyl5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Tspy15 testis-specific protein, Y-encoded-like 5 [Mus musculus (house mouse)]
Official Symbol	TSPYL5
Synonyms	TSPYL5; testis-specific protein, Y-encoded-like 5; 6330530B20; E130308C19Rik; testis-specific Y-encoded-like protein 5; TSPY-like protein 5;
Entrez Gene ID	239364
Protein Refseq	NP_001078890
UniProt ID	Q69ZB3
Chromosome Location	15 B3.1; 15
Function	molecular_function;